
COLPHARMA ®

Aneroid KIT

Sphygmomanometer
with stethoscope

12489

EN

IT

DE

FR

Fig. A / Abb. A


Fig. B / Abb. B



Fig. C / Abb. C

Fig. D / Abb. D

Aneroid KIT

EN

Aneroid Sphygmomanometer with stethoscope

- | | |
|---------------------------------|-----------------------------|
| 1. Ear Pieces | 7. Chest piece |
| 2. Stethoscope | 8. Arm Cuff |
| 3. Metal ear tube | 9. Y-Shaped adjustment tool |
| 4. Inflation bulb | 10. Manometer/Cuff Tube |
| 5. Monometer gauge | 11. Carrying Pouch |
| 6. Adjustable air release valve | 12. Air Intake Hole |

Read the instructions carefully before using this device.

This manual contains important information on the operation and safety of this device.

Please, read the instructions carefully before using the device and keep them handy for future reference.

The aneroid sphygmomanometer is a device used for measuring blood pressure through auscultation based on the Korotkoff method.

OVERVIEW

The Colpharma® aneroid sphygmomanometer features an integrated stethoscope and a self-adjusting cuff with a D-ring and artery position indicators, making it perfect for self-monitoring blood pressure.

All the air chambers, including the inflation bulb (4), are in PVC.

The inflation bulb (4) also includes a built-in air release valve (6).

The manometer's dial (5) includes an indicator displaying normal blood pressure values (numbers in green) for ease of reference.

WARNING!

- **It is forbidden to use this device for purposes beyond those outlined in this manual.**
- **Colpharma Srl shall not be held liable for any damages arising from improper use or from usage that deviates from the instructions provided in this manual.**
- **The manufacturer reserves the right to modify both the device and the manual without providing advance notice with the goal of enhancing the product's characteristics.**

KEY NOTES ON BLOOD PRESSURE INFORMATION

- Blood pressure is the force exerted by the blood within the arteries, originating from the heartbeat. Two key values are measured: the systolic (maximum) and diastolic (minimum) pressure.
- Consistently high blood pressure values require a medical consultation to prevent serious health issues.
- Regularly sharing your blood pressure readings with your doctor and bringing any anomalies to his/her attention is very important. Never base your conclusions on a single measurement.
- A significant increase in blood pressure may have various potential causes, which your doctor will investigate, explain and address appropriately.
- Never change the prescribed medication dosage or initiate medical treatments without prior medical consultation.
- Blood pressure can fluctuate significantly throughout the day, influenced by physical activity and overall conditions. Measurements should ideally be taken when resting or relaxing. At least two readings are recommended every time (once in the morning before taking any medication or food, and once in the evening before bedtime, taking a shower, or any medications) calculating the average of the measurements taken.
- It is normal for two closely spaced readings to differ.
- Readings taken at the doctor's clinic or at a pharmacy location

may vary from those taken at home because the environments are different.

- Repeated measurements provide more reliable data than single measurements.
- Wait at least 5 minutes between measurements.
- Regular blood pressure monitoring during pregnancy is very important because significant changes may occur.

CHECKING YOUR BLOOD PRESSURE

Here is a table showing adult blood pressure values according to international guidelines (ESH, ESC, and JSH). The data is expressed in mmHg. The assessment is based on the higher value, for instance, 140/80 mmHg or 130/90 mmHg indicates elevated blood pressure.

| Type | Systolic | Diastolic | Recommendations |
|------------------------------------|----------|-----------|---------------------------|
| 1. Low blood pressure | ↓ 100 | ↓ 60 | See a doctor |
| 2. Optimal blood pressure | 100-130 | 60-80 | Self-monitoring |
| 3. Elevated blood pressure | 130-135 | 80-85 | Self-monitoring |
| 4. High blood pressure | 135-160 | 85-100 | See a doctor |
| 5. Dangerously high blood pressure | ↑ 160 | ↑ 100 | See a doctor immediately! |

MEASURING BLOOD PRESSURE

KEY INSTRUCTIONS FOR A RELIABLE MEASUREMENT

1. Avoid smoking, consuming caffeinated or alcoholic beverages, eating, and exercising 30 minutes before taking the measurement.
2. Empty the bladder and stay in a quiet, comfortable room for 5 minutes.
3. Sit on a chair with a backrest and relax for 5 minutes before taking the measurement. Keep your feet flat on the ground and avoid crossing your legs.
4. Remove any bulky or thick clothing that may constrict the arm. To avoid constriction, do not roll up your shirt sleeves – they will not interfere with the cuff if worn over them.
5. Ensure that the cuff is positioned 1-2 cm above the elbow (Fig. B.1).
6. Check that the cuff is at about heart level (Fig. B.2) and aligned with the brachial artery (ARTERY) (Fig. B.3).
7. **Avoid tightening the cuff (8) excessively, as this can significantly restrict blood circulation.**
8. To ensure an accurate reading, properly regulate the deflation process. It is advisable to deflate at a rate of 2-3 mmHg per second, or a decrease of 1-2 notches on the manometer (5).

IMPORTANT: Do not press the stethoscope too firmly against the cuff, as this could lead to inaccurate diastolic (minimum) pressure measurements.

WARNING!

The reading from this device is not a diagnosis.

It cannot replace the advice of your primary care provider, particularly if the results do not match your symptoms.

It is important not to depend exclusively on this measurement; also take into account any other symptoms and the overall condition of the patient.

If needed, call a doctor or an ambulance.

MEASURING PROCEDURE

1. Position the chest piece (7) beneath the cuff (8) or slightly lower it at a distance of 1-2 cm, ensuring it is in contact with the skin and over the brachial artery.
2. Insert the ear pieces (1) into your ears, double-checking that the chest piece (7) is correctly positioned for clear hearing of the Korotkoff sounds.
3. Remain stationary during the blood pressure measuring process.
3. Secure the inflation bulb air release valve (6) by rotating its knob clockwise, being careful not to overtighten it.
4. Inflate the cuff (8) to a level about 40 mmHg above the expected systolic (maximum) pressure. When self-measuring, use the inflation bulb (4) with your free hand and avoid using the hand from the arm being measured.

IMPORTANT: If the expected blood pressure value is unknown, inflate the cuff to 200 mmHg.

5. Slowly open the air release valve (6) by turning the adjustment wheel counterclockwise to allow a consistent pressure drop of about 2-3 mmHg per second. Pay attention to notice when the cuff (8) starts deflating. Record the systolic value when you first hear faint, rhythmic, or strong sounds.
6. Continue to steadily decrease the pressure of the cuff at the same speed. When you distinctly hear at least two consecutive heartbeats, that indicates the maximum (SYSTOLIC) pressure, while when these heart sounds become impossible to hear you have found the minimum (DIASTOLIC) pressure.
7. Fully open the air release valve (6) to deflate the cuff (8).
8. Avoid keeping the cuff (8) inflated for an extended period.
9. Repeat the measurement at least twice, writing down the values, date, and time right after each measurement.
10. Remove the cuff (8) and the stethoscope (2).
11. Store the meter with the valve (6) set to the fully open position.

CALIBRATION AND ACCURACY TESTING

WARNING!

Check the accuracy of this device every two years or after any mechanical impact (such as a fall) to ensure accurate readings. The meter's accuracy can be checked visually. Just make sure the manometer's needle (5) rests on 0 mmHg in the dial when the unit is completely deflated. If the needle of the manometer (5) deviates from the 0 mmHg reference point, calibration is required (see Fig. C).

IMPORTANT: Tolerance Range: +/- 3 mmHg from the center of the 0 mmHg position. Recalibration is needed if the needle falls outside this range.

- Gently disconnect the PVC tube (10) from the manometer (5), paying attention not to damage the tube.
- Place the Y-shaped tool (9) provided to adjust the manometer's needle (5), placing it around the air intake hole (12).
- Rotate it gently clockwise to move the needle (5) left. Subsequently, rotate the tool (9) counterclockwise to shift the needle (5) to the right.
- Align the needle back to the 0 mmHg mark.
- After adjustment, remove the Y-shaped tool (9).
- Gently reconnect the tube (10) to the air intake hole (12).

WARNING!

Avoid adjusting the needle inside the manometer (5) with the Y-shaped tool (9) beyond the specified tolerance range of +/- 3 mmHg. Otherwise the device will get damaged. If the device remains unstable post-calibration, contact your retailer.

For details on coverage, consult the warranty.

TROUBLESHOOTING/RESOLVING MALFUNCTIONS

In case of operational issues, perform the following checks and adopted the recommended solutions:

Weak, altered, or interfered heartbeat sound perception:

- Check the ear pieces (1) and ear tube (3) for cleanliness and cracks. Make sure they are correctly fitted in your ears.
- Check the tube (10) for any breaks or twists.
- Ensure the chest piece (7) is undamaged and properly positioned over the skin in the brachial artery position (see Fig. B.3).

The pressure does not increase despite Inflation using the inflation bulb (4):

- Check that the air release valve (6) is closed.
- Check that the connections among the cuff (8), inflation bulb (4), and manometer (5) are correct.
- Inspect the cuff (8), tube (10), and/or inflation bulb (4) for any damage or air leaks.

Valve (6) does not allow setting the deflation speed at around 2-3 mmHg/sec.:

- Remove the air release valve (6) from the inflation bulb (4) and check for obstructions or dirt.
- Clean or remove any obstructions and retry.
If the issue persists, contact your retailer.

The manometer needle (5) is not at 0 ± 3 mmHg when in the rest position:

- Ensure the air release valve (6) on the inflation bulb (4) is fully open.
- If the deviation is more than 3 mmHg, recalibrate the manometer using the Y-shaped tool (9) provided (refer to the “CALIBRATION AND ACCURACY TESTING” section in this manual). If the issue persists, contact your retailer.

CARE AND MAINTENANCE

The Colpharma® aneroid sphygmomanometer is **EC** approved and provided pre-calibrated and thoroughly inspected before shipping.

Prior to each use, especially when at rest, check that the aneroid sphygmomanometer's needle on the dial is within the specified tolerance range of +/- 3mmHg.

Keep the device in its provided pouch (11) to keep it clean and protected.

CLEANING

MANOMETER - Only use a soft, dry cloth for cleaning.

INFLATION BULB - Wipe with a slightly damp and soft cloth and then dry it completely.

CUFF - To clean the cuff properly, gently remove any dirt using a cloth dampened with soap and water.

STETHOSCOPE - Clean with a moistened, soft cloth and then let it dry completely. For other components (ear pieces and chest piece) clean using lukewarm soapy water dry completely.

SAFETY

- This device must be used exclusively as described in this manual. The manufacturer shall not be liable for any damages caused by improper use.
- Never disassemble the manometer (5).
- This device has been manufactured using sensitive components and it must be handled with care. Please, follow the storage and operation instructions detailed in the “TECHNICAL SPECIFICATIONS” chapter.
- Protect the device from strong shocks and impacts. Shocks and vibrations may alter the device calibration, leading to inaccurate readings.
- Avoid damaging the PVC components using sharp or pointed items.
- Never prod the cuff fabric or parts (8) with sharp objects. Doing so may damage the cuff and parts.
- Do not store the device in high-humidity environments
- Do not drop the device.
- Protect the device from direct sunlight.
- Keep the device away from heat sources.
- Do not use the device if you deem it damaged or if you notice any anomaly.
- Only inflate the cuff after wearing it.
- Do not exceed 300 mmHg when inflating.
- Deflate the cuff completely after every use.
- Measurement time should not exceed 2 minutes.
- Read the additional safety information in the various sections of this manual.

WARNING!

Children must always use this device under adult supervision. Some parts are small and may represent a choking hazard if swallowed.

Cables and tubes may represent a choking hazard.

RECOMMENDATIONS

- Avoid using this device in moving vehicles (such as cars or planes).
- Do not use this device with cuffs for newborns.
- Knowledge and regular training are key to prevent reading errors.
- Be aware that noisy environments may disrupt measurement accuracy.
- Skilful handling is required to ensure the cuff deflates at the correct speed.
- Accurate measurements require excellent eyesight and hearing.

DISPOSAL GUIDELINES

Do not dispose of this device using standard municipal solid waste services. Responsibly dispose of the device at local community recycling centers, ensuring proper recycling of its materials.

SYMBOL

Lot number

Read the instructions

Caution - Safety Instructions

Compliant with European regulation

Medical device

Manufacturer

Imported and distributed by

Authorised Representative for Switzerland

Authorised Representative for the European Community

WARRANTY

This device is covered by a 2-year warranty from the date of purchase, against manufacturing defects, valid only upon submission of the purchase receipt.

The warranty does not include:

- costs and risks associated with transportation.
- damage due to improper use or failure to follow the operating instructions.
- damage from falls or misuse.
- packaging/storage materials and operating instructions.
- regular inspections and maintenance (calibration).

Note that repairs or replacements under warranty do not extend or renew the original warranty period.

TECHNICAL SPECIFICATIONS

| | |
|-----------------------|---|
| Dial diameter | 48 mm |
| Manometer dimensions: | 50 * 50 * 24 mm |
| Product net weight: | 105 ± 5 g |
| Measurement range: | 0-300 mmHg |
| Resolution: | 2 mmHg |
| Accuracy: | ±3 mmHg |
| Air leakage: | ≤4 mmHg/min |
| Deflation speed: | 2-3 mmHg/s |
| Rapid deflation: | from 260 mmHg to 15 mmHg ≤ 10s |
| Operating conditions: | temperature +10 °C - +40 °C; humidity 20% - 85% |
| Storage conditions: | temperature -20 °C - +70 °C; humidity: 20% - 85% |

We reserve the right to apply technical changes.

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HONSUN (NANTONG) Co., Ltd
Shanghai International Holding Corp.
GmbH (Europe) Eiffestrasse 80,
MedNet SWISS GmbH,